



www.bunsen.es



SAND BATHS

Ref. 1560/1561

GENERAL CHARACTERISTICS

Equipment designed to perform, in dry, tests at high temperatures.

- * Stainless steel surface, easy to keep clean and resistant to corrosion AISI 316.
- * Heating elements distributed evenly throughout the whole plate thus providing good uniformity in temperature.
- * AISI 304 stainless steel basin, with easy positioning and whose absence allows the use of the bath as a heating plate.
 - * Metallic cabinet painted in epoxy.
- * Excellent thermal isolation that prevents the transmission of heat to the cabinet where the control elements are located.
- * TERM-OFF safety system. Which limits the maximum temperature of the equipment, thus prolonging the useful life of the same.
 - * Maximum temperature on the surface of the plate up to 400°C.
 - * Electronic Voltage regulation by means of energy pulses.

Manufactured under EC directives.



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN

Tlf. +34 916 113 584 - Fax. +34 916 128 254

info@bunsen.es - www.bunsen.es



| Ref | Useful measures W x H x D (mm) | Exterior measures W x H x D (mm) | Power (W) |
|------|--------------------------------|----------------------------------|-----------|
| 1560 | 200 x 60 x 400 | 190 x 220 x 480 | 2000 |
| 1561 | 250 x 60 x 500 | 190 x 280 x 580 | 3000 |

| Ref | Accessories |
|------|---|
| 1735 | Digital programmer |
| 1521 | Power outlet to connect a contact thermometer or a digital programmer |



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
 Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es

