



www.bunsen.es



DISTILLER 4L.

Ref. DA1040

GENERAL CHARACTERISTICS

- * Automatic distiller for continuous production. Excellent quality of distillation.
- * Furniture, cooling and heating elements made of stainless steel. Housing made of stainless steel, painted with epoxy.
 - * Easily removable for cleaning
 - * With thermostat, that protects the distiller in case of overheating.
- * Safety mechanism which switches of the heaters in the absence of cooling water and connects it again.
- * It includes a system for limiting the flow of water from the tap water, optimizing the flow down to a necessary minimum.
 - * Adaptable water inlet connection to the inlet pipe.
 - * Output connector adaptable to a vessel.
 - * Water consumption; 1L/min.
 - * Conductivity at 20°C ($\mu\text{S/cm}$); 1,5.

Manufactured according to EC directives.

Ref	Capacity	Voltage	Dimensions W x H x D (mm)	Weight	Power (W)
DA1040	4 L/h	230 V. 50/60 Hz. 13 A	370 x 435 x 220	8 Kg.	3000



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

