



www.bunsen.es



DISTILLER WITH TANK

Ref. DD1048

GENERAL CHARACTERISTICS

- * Automatic distiller and continuous output.
- * Built-in water storage tank with a capacity of 16 litres.
- * Cooling element and heaters made of stainless steel.
- * Excellent quality of distillation. Dispenser tap of the distilled water.
- * Safety device protects the water still in the event of water supply failure.
- * Some models have a device to fill a container without worrying when it finished. It cuts water and electrical supply to avoid useless expenses of water and electricity.
 - * With thermostat, that protects the distiller in case of overheating.
- * Flow restrictor that limits the water flow from tap water, optimizing it to a necessary minimum.
 - * Water consumption; 0,9 L/min.
 - * Conductivity at 20°C (µS/cm); 1,5.

Manufactured according to EC directives.

Ref	Capacity	Voltage	Dimensions W x H x D (mm)	Tank	Weight	Power (W)
DD1048	8 L/h	230 V. 50/60 Hz. 26 A	710 x 593 x 340	16L.	34 Kg.	6000



BUNSEN SA. C/ Tokio, 2 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

