



www.bunsen.es



BATH FOR OIL SAMPLING 4C.

Ref. 1625-4

GENERAL CHARACTERISTICS

- * Stainless steel box AISI 304.
- * Aluminium block fixed.
- * Heating by means of a resistance that cover all surface of the aluminium block.
- * Digital temperature control by digital thermostat controlled by microcontroller from environment + 5° to 30°C. Resolution, 0.1°C.
- * 3 digits temperature display.
- * Specially cavities designed by International Olive Council (IOC) according to EU N° L.248/67 standards.

Manufactured under EC directives.



Ref	Probe	Power	Block dimensions W x D	Exterior dimensions W x H x D	Max. temp.	Precision °C
1625-4	PTC	75W.	150 x 150mm.	185 x 80 x 230mm.	30°C	±0.5%



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

