



[www.bunsen.es](http://www.bunsen.es)



# BATH FOR OIL SAMPLING 6C.

Ref. 1625-6

## GENERAL CHARACTERISTICS

- \* Stainless steel box AISI 304.
- \* Aluminium block fixed.
- \* Heating by means of a resistance that cover all surface of the aluminium block.
- \* Digital temperature control by digital thermostat controlled by microcontroller from environment + 5° to 30°C. Resolution, 0.1°C.
- \* 3 digits temperature display.
- \* Specially cavities designed by International Olive Council (IOC) according to EU N° L.248/67 standars.

Manufactured under EC directives.

Más baños



Other baths

Ref	Probe	Power	Block dimensions W x D	Exterior dimensions W x H x D	Max. temp.	Precision °C
1625-6	PTC	125W.	150 x 230mm.	185 x 80 x 315mm.	30°C	±0.5%



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fa. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

