



www.bunsen.es



AUTOMATIC DISTILLER 4L.

Ref. DA1148

GENERAL CHARACTERISTICS

- * Automatic distiller and continuous output..
 - * High quality distillation with water of high purity and free of metal ions and pyrogens.
 - * Safety device protects the water still in the event of water supply failure.
 - * The different parts made of borosilicate glass and heater protected by a quartz sheath lengthens its duration.
 - * Small size for easy handling and cleaning. All models feature with wall mounting bracket.
 - * Water consumption; 1 L/min.
- Conductivity at 20°C ($\mu\text{S/cm}$); 1.

Manufactured according to EC directives.

Ref	Capacity	Voltage	Dimensions W x H x D (mm)	Weight	Power (W)
DA1148	4.2 L/h	230 V. 50/60 Hz. 13 A	590 x 510 x 170	5.2 Kg.	3000



BUNSEN SA. C/ Tokio, 2 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

