



[www.bunsen.es](http://www.bunsen.es)



# DISTILLER 8L.

Ref. DA1048

## GENERAL CHARACTERISTICS

- \* Automatic distiller for continuous production. Excellent quality of distillation.
- \* Furniture, cooling and heating elements made of stainless steel. Housing made of stainless steel, painted with epoxy.
  - \* Easily removable for cleaning.
- \* With thermostat, that protects the distiller in case of overheating.
- \* Safety mechanism which switches of the heaters in the absence of cooling water and connects it again.
- \* It includes a system for limiting the flow of water from the tap water, optimizing the flow down to a necessary minimum.
  - \* Adaptable water inlet connection to the inlet pipe.
  - \* Output connector adaptable to a vessel.
  - \* Water consumption; 1,4 L/min.
  - \* Conductivity at 20°C (μS/cm); 1,5

Manufactured according to EC directives.

Ref	Capacity	Voltage	Dimensions W x H x D (mm)	Weight	Power (W)
DA1048	8 L/h	230 V. 50/60 Hz. 26 A	370 x 635 x 260	14 Kg.	6000



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fa. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

